

medbar[®]

Pouch Perforator



REF 286 02



CE 2292

STERILE EO



Pouch Perforator

Intended Use

- Safe, disposable perforator used to rupture the maternal amniotic sac in order to induce labor, i.e., perform Amniotomy (a mechanical method for the induction of labor).

Product Features

- The specially designed tip allows the amniotic membrane(pouch) to be perforated in a single motion.
- Rounded tip designed to rupture the amniotic membrane without injury to the mother or fetus.
- Total length is 26 cm.

Quality

- Manufactured under ISO 13485:2012 Quality Management Standard.
- CE marked according to 93/42/EEC Medical Device Directive and classified as Class I sterile medical device.
- Notified Body: UDEM Ltd. Şti. (2292)
- EC Certification No: M.2016.106.7000.
- GMDN Code of Insemination Cannulas: 46520.

Sterilization

- Insemination Cannulas sterilized with Ethylene Oxide.
- Validated according to ISO 11135 standard.

Biocompatibility

- Latex-free.

Shelf Life

- 5 years.



REF NO	PRODUCT NAME	STERILIZATION
210 01	Pouch Perforator	EO

